THE EVOLUTION OF ANL CMT GLOVEBOXES

by

Richard F. Malecha, Arthur A. Frigo, and Daniel E. Preuss

Chemical Technology Division Argonne National Laboratory 9700 South Cass Avenue Argonne, Illinois 60439

Telephone: (630)252-4510 Fax: (630)252-5735

Submitted for presentation at:

2000 Annual Conference and Equipment Exhibit of the American Glovebox Society August 21-23, 2000 New Orleans, Louisiana

The submitted manuscript has been created by the University of Chicago as Operator of Argonne National Laboratory ("Argonne") under Contract No. W-31-109-ENG-38 with the U.S. Department of Energy. The U.S. Government retains for itself, and others acting on its behalf, a paid-up, nonexclusive, irrevocable worldwide license in said article to reproduce, prepare derivative works, distribute copies to the public, and perform publicly and display publicly, by or on behalf of the Government.

Work supported by the U.S. Department of Energy under Contract W-31-109-ENG-38.

INTRODUCTION

- Argonne National Laboratory
- Chemical Technology Division (CMT)
- Modular Gloveboxes Designed for Experimental Work
 - Laboratory Scale
 - Engineering Scale
- Experimental Work Includes Equipment and Process Development for:
 - Pyrochemical
 - Nuclear Waste Treatment

Argonne National Laboratory

Electrochemistry



HISTORY

- The First ANL-CMT Division Modular Glovebox was Designed and Installed in 1959. It was designated as the CENHAM glovebox.
- Design Objective:
 - Provide Modular Configuration
 - Provide Controlled Atmosphere Environment
 - Maximize Viewing Accessibility
 - Provide "Standardized" Work Area for Laboratory Research Work
 - Provide Modular Utility Service Access
 - Include User-Friendly Considerations



GLOVEBOX DESIGN REQUIRES AWARENESS OF USER NEEDS

- Evolution of Design Over the Past 40 Years
- Request for New Gloveboxes Usually Based Upon Similar Attributes of an Existing Enclosure Plus Particular Project Changes for Use.
- Discussion with Staff and Laboratory Operating Personnel.
- Considerations:
 - Simplicity in Design
 - Cost Saving
 - Schedule
- Defining "What is Necessary" and "What Would be Nice."

Argonne National Laboratory



USER-FRIENDLY GLOVEBOX DESIGN

- Gloveboxes Use a Modular Design Concept
- Glovebox Size is Designated as Modules in Length and Tiers in Height
- Basic Module is 42 inch Cube
- Modular End Plates Provide for:
 - Utility Services
 - Ventilation/Purification System

Argonne National Laboratory

- Filter Housings
- Transfer Locks
- Bagports



MATERIALS OF CONSTRUCTION

- Steel Shell (Painted)
- Stainless Steel Shell
- Steel Support Frame

Argonne National Laboratory

- Gloveports
- Glass Windows



Argonne National Laboratory

STRUCTURE

- Glovebox Shell
- Structural Elements
- Unistrut
- Floor
- End Plates
- Hoists
- Floor Wells



UTILITIES

Process Feedthroughs for Services

Argonne National Laboratory

- Electrical
- Instrument
- Gas or Liquid
- Lighting



TRANSFER SYSTEMS

- Bagout
- Large Horizontal Transfer Lock
- Small Horizontal Transfer Lock

Argonne National Laboratory

- Vertical Transfer Lock
- Sphincter



QA APPLICATION

- Management Control
 - Team Approach
 - Process Development
 - Design and Fabrication
 - Scheduling
- Design Review
- Design for Functionality and Manufacturability Reviews (DFM)
- System Design Description
- Safety Review
- Operational Readiness Review



WINDOWS

- Window Viewing Area Comprises Approximately 60% of the Glovebox Side Walls that Suffices for Monitoring Experimental Equipment and Process Operation
- Weatherstrip Type "Zipper" Seal Used for Window Installation
- Window Concept
 - Nominal 36 Inch Square Windows and Window Openings with Rounded Corners
 - 3/8 Inch Thick Laminated Safety Glass Windows
 - Gloveports are Attached Through the Windows







RFMAGS2000_12





RFMAGS2000_13

GLOVEBOX SUPPORT STAND





RFMAGS2000_14

4 MODULE - 1 1/2 TIER GLOVEBOX





RFMAGS2000_15

THE EVOLUTION OF ANL CMT GLOVEBOXES







RFMAGS2000_16

GAS RECIRCULATION SYSTEM





RFMAGS2000_17





RFMAGS2000_18

GLOVEBOX FILTER ACCESSIBILITY





RFMAGS2000_19

EPOXY RESIN

- Inside Weld Joints are Caulked with an Epoxy Resin.
- The Epoxy Resin (ABAWELD) has been Used on Gloveboxes in CMT for Many Years
- The Caulked Joints Provide a Smooth Corner Fillet Designed for Easy Clean-Up

Argonne National Laboratory



WEATHERSTRIP INSTALLATION





RFMAGS2000_21

GLASS TYPE

- Glass laminate per ASTM C1172-91
- Two lite laminate of Kind LA, Class 1, q³ quality glass
- Each lite is .19 inch thick with overall composite thickness of .38 inch

Argonne National Laboratory



FAILURE OF GLOVEBOX WINDOW WITH GLOVEPORTS





RFMAGS2000_23

FAILURE OF GLOVEBOX WINDOW WITHOUT GLOVEPORTS





WINDOW TEST RESULTS

- Window with four gloveports
 - Test terminated due to structural failure of the window glass
 - Cracks developed in the glass extending across the two lower gloveports
 - Test pressure at failure was 17 inches of water
 - No evidence of seal failure

Argonne National Laboratory



WINDOW TEST RESULTS

- Window without gloveports
 - Test terminated due to structural failure of the window glass
 - Crack developed in the glass radiating out from the center
 - Test pressure at failure was 30 inches of water (> 2 psig)
 - No evidence of seal failure

Argonne National Laboratory



INTERIOR ENVIRONMENT ATMOSPHERE

- Flow Controls Once Through Gas
- Inert Gas Recirculation System -Requires Purification System with Filters

Argonne National Laboratory

Pressure Controls



SUMMARY

- Design Approach Based Upon User-Friendly Concept
- Utilization of Existing Component Designs
- Cost Effective
- Schedule
- Adaptable to Project Process Changes Without Losing Overall Effectiveness of "User-Friendly" Approach.

